

ATMOSPHERIC CHEMISTRY RESEARCH AT THE AIR POLLUTION RESEARCH CENTER, UCR

Roger Atkinson, Janet Arey, Jingsong Zhang
and Paul Ziemann

CDAWG

Big Bear, October 20, 2006

Photochemical Air Pollution

- Photochemical air pollution is caused by the interaction of sunlight with oxides of nitrogen ($\text{NO}_x = \text{NO} + \text{NO}_2$) and volatile organic compounds (VOCs) emitted into the atmosphere.

Emissions of VOCs and NO_x

- ~1,500 Million tons per year of organic compounds (not including methane) are emitted into the atmosphere from natural (or biogenic) sources and from human activities.
- On a global basis, 80-90% of these non-methane organic emissions are from vegetation.
- NO_x is emitted mainly from combustion sources – vehicles and fossil-fueled power plants being prime examples.

Biogenic VOCs from vegetation



Anthropogenic emissions



Consequences of VOC emissions

In the presence of oxides of nitrogen (NO_x) and sunlight, a series of reactions occur leading to photochemical air pollution, which includes:

- Formation of ozone and other toxic chemicals.
- Formation of secondary aerosol (in addition to directly emitted aerosol).




Morning, vehicle
emissions

Afternoon,
pollution



50 km altitude




transported downward
by convection



STRATOSPHERE

15 km



TROPOSPHERE

The hydroxyl (OH) radical is a key reactive intermediate

The atmosphere is a large chemical reactor




Current research at APRC

- Kinetics, products and mechanisms of gas-phase reactions of individual VOCs with OH radicals, NO₃ radicals and O₃.
- Chemical pathways leading to formation of secondary organic aerosol from reactions of individual VOCs.
- Ambient measurements of selected organics, including PAHs and nitro-PAHs, and comparison with laboratory studies.
- Real-time ambient measurements of selected chemical species by spectroscopic methods.

Gas-phase studies

Janet Arey


Roger Atkinson



One of three ~7000 liter
volume Teflon chambers.

20% of blacklamps are on.


5870 liter volume evacuable chamber



5870 liter volume evacuable chamber






Direct air sampling atmospheric pressure ionization tandem mass spectrometer (API-MS/MS)




Gas and aerosol studies: alkanes

- Liquified petroleum gas (LPG).
- 50-60% of gasoline.
- 50-60% of organics in vehicle exhaust.
- ~50% of non-methane VOCs in air in urban areas.



Alkanes react with OH radicals in the atmosphere




Product yields (as of 1982-1994)

alkane	% carbonyls	% nitrates	% missing
<i>n</i> -butane	84	8	8
<i>n</i> -pentane	~30	12	~58
<i>n</i> -hexane	low	21	high

From data in Carter *et al.*, 1976; Atkinson *et al.*, 1982.


Missing products postulated to be 1,4-hydroxycarbonyls, formed after alkoxy radical isomerization (Carter *et al.*, 1976)



Further progress was halted until new analytical techniques became available




API-MS analyses in positive ion mode



API-MS spectra of
OH + *n*-butane, *n*-
pentane and *n*-
hexane


API-MS analyses in negative ion mode



OH + C₅-C₈
n-alkanes:
NO₂⁻ adducts


Solid-Phase MicroExtraction fibers

Can be coated with a chemical for on-fiber derivatization of carbonyl-containing compounds, followed by GC analyses.



Fibers are retractable. After coating with PFBHA, the fibers are exposed to the chamber contents with the chamber mixing fan on for typically 5 min.

Fiber, ~20 mm long



Product yields (%) in air at 50% RH

alkane	carbonyls	nitrates	hydroxy-carbonyls	hydroxy-nitrates
<i>n</i> -pentane	50	10.5	54 (36)	2.6
<i>n</i> -hexane	10	14	57 (53)	4.6
<i>n</i> -heptane	<1	18	51 (46)	4.7
<i>n</i> -octane	<1	23	53 (27)	5.4


Yields in parentheses from API-MS in dry air (3-5% RH)

Success so far!



But what happens to 1,4-hydroxycarbonyls in the atmosphere?


5-Hydroxy-2-pentanone (5H2PO) is the only commercially available 1,4-hydroxycarbonyl



FT-IR spectra of a
5H2PO - N₂
mixture in the dark

5-Hydroxy-2-pentanone has lifetime of 1.1 hr in dry air.

4,5Dihydro-2-methylfuran is stable in dry air, but converts to 5H2PO at 5% relative humidity, with a lifetime of 3.5 hr.




Reactions of 5-hydroxy-2-pentanone and 4,5-dihydro-2-methylfuran

- 4,5DH2MF reacts very rapidly with OH radicals, NO₃ radicals, and O₃. Daytime lifetimes:

<i>n</i> -Pentane:	1.5 day	(OH)
5-Hydroxy-2-pentanone:	9 hr	(OH)
4,5-Dihydro-2-methylfuran:	38 min	(OH)
	7 min	(O ₃)


Relative humidities where dihydrofuran formation is important





Aerosol-phase studies

Paul Ziemann




Experimental Apparatus



Real-Time Mass Spectra of SOA from Decane and Pentadecane + OH/NO_x



Reactions of 1,4-Hydroxycarbonyls



substituted
tetrahydrofuran



Atmospheric lifetimes

$$\tau_{OH} \sim 1.3 \text{ h}$$




$$\tau_{NO_3} \sim 24 \text{ s}$$

$$\tau_{O_3} \sim 7 \text{ min}$$

Desorption Profiles of SOA from Decane and Pentadecane + OH/NO_x




SOA Yields from Linear Alkanes + OH/NO_x




Zingsong Zhang

Spectroscopic studies


Cavity Ring-Down Spectroscopy (CRDS)



- Quantitative and absolute measurements
- High sensitivity with a long effective absorption path ($\sim \text{km}$) in a compact setup ($\sim \text{m}$)
- Spectroscopic selectivity
- Real-time and in-situ detection




Typical CRDS Signal



Usable data in $\sim 3.5 \tau$, effective absorption length $\sim 21 \text{ km}$ (Speed of light = $300 \text{ m}/\mu\text{s}$)


Cavity Ringdown Spectrometer



Mirror


Mirror

CRDS Absorption Spectrum of HONO



HONO Concentration ~ 1 ppm


CRDS Detection of HONO



Linear response
Large dynamic range
High detection sensitivity


Number density of HONO is based on measurements using a NO-NO_x Analyzer
1 ppm = 2.38×10^{13} molecules/ cm^3 at 735 Torr and 25°C

Detection Sensitivity for HONO




- Detection sensitivity **5 ppb/15 s (3σ)**
- Recently improved to **1 ppb/15 s (3σ)** with high-quality mirrors ($\tau_0 \sim 10 \mu\text{s}$)
- Ambient HONO **0.1 ppb – 15 ppb**


CRDS Spectrum of NO₂



CRDS Detection Sensitivity for NO₂



Variations in NO₂ concentration during daytime





Research funded by:

California Air Resources Board

U.S. Environmental Protection
Agency

National Science Foundation

U.S. Department of Energy

Thank you for your attention

